

## 2 Arms LED Illuminator



This is a 6 Watt LED dual light microscope illuminator in mint condition, with two long LED cable that can be bent into various positions. The light source is: 100-240VAC / 50-60Hz. Unit. The lamps are 3Watt each (6 Watt in total) with fully adjustable intensity.

Perfect for Evenly Illuminating a Specimen without Heating This high intensity cool lighting LED system is excellent for all applications including biological and electronics applications.

Configurable LED Peak Wavelength or Color Temperature From UV, Visible to NIR or Cool / Neutral / Warm white

6400K white color

Lamp Color Temperature: 3250 Degrees Kelvin

LED power: 3WX2pcs

Input Voltage: 85-256VAC / 50-60Hz

Fuse: 1A / 250V

60 bright white LED bulbs

100,000 hours of life

Intense and focused shadow-free illumination

The brightness of 0-100% Adjustable, high brightness

Suit for the dia of the lens: 23mm-61mm

Power: 3W

Dimension:140mm x 106mm x 70mm,(LxWxH )

Long dual gooseneck lights:590mm/each

colour:black

Material:metal

Package Contents:

2 \* Gooseneck LED

1 \* Adapter EU

1 \* Controller

Plug Type: EU Plug, it applies to Indonesia, Germany, France, Spain, Italy and so on.

Only the above package content, other products are not included.