


## MPlanFL 50x/0.80 BD DIC $\infty$ /0 WD: 1.0

MPlanFL 50x/0.80 BD DIC  $\infty$ /0 WD: 1.0 is BD Semi-Apochromatic Objective: A lab-grade all rounder delivering high resolution. Combines dual-mode imaging with semi-apochromatic correction, capturing intricate material details while ensuring faithful color reproduction-meeting the demands of advanced research and high-end quality inspections.



Objective	Magnification	NA (mm)	Working Distance (mm)	Field of View	Barrel
MPlanFL 50x/0.80 BD DIC $\infty$ /0 WD: 1.0	50X	0.80	1.00	25	180